

How to Collect a Pinworm Sample

Read all the instructions before collecting a sample. It's important to follow the instructions so the sample is collected the right way. This will help the test results be correct.

How to Collect a Sample

- Collect the sample in the morning before bathing or using the toilet.
- Take the clear plastic paddle attached to the cap out of the collection container. Do not touch the paddle.
- The sticky surface of the paddle will be marked "sticky side."

Cómo tomar una muestra de lombrices intestinales

Lea todas las instrucciones antes de tomar una muestra. Es importante seguir las instrucciones para que la muestra se tome de la manera correcta. Esto ayudará a que los resultados de la prueba sean correctos.

Cómo tomar una muestra

- Tome la muestra por la mañana antes de bañarse o ir al baño.
- Retire la paleta de plástico transparente que se adjunta a la tapa del recipiente de recolección. No toque la paleta.
- La superficie pegajosa de la paleta estará marcada como "lado pegajoso" ("sticky side.")

(over)

- Separate the buttocks. Press the sticky side of the paddle against several areas of the skin around the anus. Do not put the paddle inside the rectum.
- Put the paddle back into the container.
- Put the cap back on the container. Make sure it is closed tightly.
- Write the patient's name, date of birth and the time the sample was taken on the container label.
- Bring the container back to the clinic right away. Keep the sample at room temperature.
- If you cannot bring the container to the clinic within 24 hours, refrigerate it until you can do so.

- Separe las nalgas. Presione el lado adhesivo de la paleta contra varias áreas de la piel alrededor del ano. No ponga la paleta dentro del recto.
- Vuelva a colocar la paleta dentro del recipiente.
- Vuelva a colocar la tapa en el recipiente. Asegúrese de que esté bien cerrado.
- Escriba el nombre del / de la paciente, la fecha de nacimiento y la hora en que se tomó la muestra en la etiqueta del recipiente.
- Lleve el recipiente a la clínica de inmediato. Mantenga la muestra a temperatura ambiente.
- Si no puede llevar el recipiente a la clínica dentro de las 24 horas, refrigérelo hasta que pueda hacerlo.